




High level instructions for easy and consistent
eye dominance measurement with the
EyeDominanceApp




Have your subject stand about 8-10 feet away from you while you're holding your phone facing them.




Have them face you directly with their right hand centered on chest with the index finger extended. They should be looking directly at the camera lens on the back of the phone facing them.



Have them look directly at the camera lens and extend their hand and point without looking back at their finger




Take a photo with the right hand pointing



Take a photo with the right hand pointing



**Repeat steps for
left arm/hand.**



After capturing (or loading from photo library) pointing photos from each hand then

- Select pupils and pointing finger location in each image within app.
- Remember the “right pupil” will be on the left of the image.
- Click analyze button
- Review results (current and historical)